

DETERGENTS

VIRUDET® N

A NEUTRAL LOW-FOAM CLEANER FOR WASHER DISINFECTORS & MANUAL APPLICATIONS

Virudet® N is for use in washer/disinfectors, AER (Automated Endoscope Reprocessor), other automated decontamination systems and manual cleaning. It is a neutral, concentrated detergent with a wide material compatibility which gently removes blood, proteins, mucus and tissue residues.

This chemistry has been specially developed for the gentle cleaning of delicate surgical instruments, utensils, anaesthesia and endoscope material, allowing simultaneous reprocessing of instruments made of different materials. It has optimal cleaning performance at temperatures of 55-60°C but for endoscopes temperature should not exceed 60°C and can be used for the following: stainless steel instruments (scalpels, tweezers etc.), articulated instruments (scissors, clamps etc.), ward utensils (bowls, feeding bottles etc.), anaesthesia utensils (breathing balloons, tubes etc.), shoes of the surgical team, transport containers made of stainless steel as well as aluminium and endoscopes.



Primary features attributable to Virudet® N are:

- Neutral low-foam detergent cleaner.
- Suitable for removal of contaminants from sensitive instrument materials.
- Uses a blend of complexing and sequestering agents, detergents and corrosion inhibitors that are effective in hard and soft water areas.
- Reduced instrument damage and repair.
- Excellent material compatibility
- Non-fixative.
- Prevents re-deposition of contaminants and water mark staining.
- Non-corrosive, will not damage delicate coatings.
- Cost effective.
- Can be used with tap water (saves DI/RO water cost).
- Can be used without neutralizer.
- Good cleaning power even at low dosage.

Packaging

- 1L container
- 5 L container
- 10 L container
- 25 L container
- 200 L drum
- 1,000 L IBC

Virudet® N is for use in washer/disinfectors, AER (Automated Endoscope Reprocessor), other automated decontamination systems and manual cleaning.

DETERGENTS

Application

Virudet® N is supplied as a concentrated product and should be diluted prior to use.

For Manual & Ultrasonic Cleaning of instruments and equipment:

For cleaning applications, dilute 1 part Virudet® N in 100 - 200 parts water (i.e. a 0.5 - 1% solution; 10ml/L).

Soak parts for between 5 and 15 minutes in the diluted Virudet® N solution at between 20°C to 40°C (allow longer times up to 15 minutes when dealing with dried on proteinaceous deposits).

Note: soak times may vary according to the temperature used and the nature of the contaminants.

After soaking, rinse off thoroughly with clean water (temperature > 40°C) or transfer to the next cleaning stage e.g. manual cleaning, washer/disinfecter or ultrasonic equipment.

For Automated Cleaning of instruments and equipment:

Injectors should be adjusted to deliver at 5 - 10ml per litre of water dependant on the expected soil loading of the contaminated equipment and the temperature of the water.

Phase	Water Type	Time	Temp	Solution	Dosage
Pre-Wash	Cold Tap Water	3 min			
Wash Step	Hot Tap Water	10 min	55-60°C	Virudet N	5-10ml/L or 0.5 - 1%
Rinse	Hot Water	1 min			
Rinse	DI or RO Water	1 min			
Thermal Disinfection	DI or RO Water	1 min	93°C	Virudet NDA	0.1 - 0.5ml/L or 0.01 - 0.05%
Drying		15 min	110-115°C		

AER

Phase	Water Type	Time	Temp	Solution	Dosage
Pre-Wash	Cold Tap Water	3 min			
Wash Step	Hot Tap Water	10 min	25-50°C	Virudet N	5-10ml/L or 0.5 - 1%
Rinse	Hot Water	2 min			
High Level Disinfection	DI or RO Water	10 min	25°C	Virusolve+ EDS Auto Amity PAA 15 Amity PAA 5 System	25ml/L (2.5%)
		5 min	30°C		10ml/L (1%) 20ml/L (2%) parts A & B
Rinse	DI or RO Water	2 min			
Rinse	DI or RO Water	2 min			
Rinse	DI or RO Water	2 min			

VIRUDET® N
A NEUTRAL LOW-FOAM CLEANER FOR
WASHER DISINFECTOR

Virudet N Properties

- Virudet® N is biodegradable.
- The product is classed as non-hazardous in diluted forms, however contact with skin and eyes should be avoided.
- See Material Safety Data Sheets for details.

Appearance	Pale brown liquid
Odour	Odourless
pH	9.8 units (neat solution)
Density g/cm³ @ 20°C	1.03 (neat)
Solubility in water	Completely miscible
Flashpoint (Abel closed cup)	None



Amity Limited, Libra House, West Street,
Worsbrough Dale, Barnsley, South Yorkshire,
S70 5PG, UK

T: +44 (0) 1226 770787

F: +44 (0) 1226 770757

E: sales@amityinternational.com

WWW.AIHC.CO.UK

© Amity Limited all rights reserved 2021

DISTRIBUTED BY: